

# LACHAT $\text{NH}_4^+$ and $\text{NO}_3^-$ REAGENTS AND STANDARDS FOR SOIL EXTRACTS

D.J. Herman

## I. $\text{NH}_4^+$ Reagents

- A. Sodium Salicylate - Nitroprusside  
Dissolve 300 g sodium salicylate and 2.0 g sodium nitroprusside in 1500 mL of deionized water, then dilute to 2 L.
- B. Sodium Hypochlorite,  
Dilute 60 mL of bleach (Approximately 5%  $\text{NaOCl}$ ) to 1 L with deionized water.
- C. EDTA Solution (60%).  
Dissolve 165g  $\text{Na}_2\text{EDTA}\cdot 2\text{H}_2\text{O}$  in about 2 L deionized water. Add  $\text{NaOH}$  pellets (about 60 of them) to pH 7.0 (the solution should be clear at that point). Dilute to 2.5 L.
- D. Buffer.  
Dissolve 100g  $\text{Na}_2\text{HPO}_4\cdot 7\text{H}_2\text{O}$  and 56 g  $\text{NaOH}$  in about 1500 mL deionized water. Dilute to 2 L.
- F. Carrier.  
Dissolve 300 g  $\text{KCl}$  in 1.7 L  $\text{H}_2\text{O}$ . Bring to 2 L volume.

## II. $\text{NO}_3^-$ Reagents

- F. Ammonium Chloride.  
Dissolve 212 g of ammonium chloride and 2.5 g of disodium ethylene diamine tetraacetate in distilled water and dilute to 2.5 L., Adjust to pH 8.5 with concentrated ammonium hydroxide. **Alternatively** for high-iron extracts: Dissolve 6.8 g imidazole in 900 mL. Bring to pH 7.5 using concentrated  $\text{HCl}$ . Add 1 mL of 2% w/v  $\text{CuSO}_4\cdot 5\text{H}_2\text{O}$  and dilute to 1 L.
- G. Color Reagent.  
To approximately 1400 mL of deionized water add 200 mL of concentrated phosphoric acid. Add 80 g sulfanilamide and dissolve completely. Dissolve 2.0 g N-1-naphthylethylenediamine dihydrochloride and dilute to 2 L. Store in dark bottle at  $12^\circ$ . Stable for 1 month.

## III. Standards

- A. Stock Standard: 50 mg  $\text{NH}_4\text{-N}$  / L plus 50 mg  $\text{NO}_3\text{-N}$  / L  
Dissolve 0.1909 g of anhydrous ammonium chloride (dried at  $105^\circ\text{C}$  2 h) and 0.36107 g potassium nitrate in 800 mL of carrier. Dilute to 1000 mL. Store in  $5^\circ$  cold room.
- B. Working Standards, Set of 7:

Pipette the indicated amount of stock standard into a 50 mL volumetric flask and dilute to volume with carrier.

Standard	$\mu\text{g N mL}^{-1}$	mL Stock
1	3.000	3.000
2	1.000	1.000
3	0.300	0.300
4	0.100	0.100
5	0.030	0.030
6	0.010	0.010
7	0.005	0.005